



# start <sup>13</sup>C Acquisition *after* <sup>1</sup>H

# in TopSpin



File....  
...new  
  
**EXPNO: 2**  
  
**Experiment:**  
  
**AA\_\_13C**



keyboard entry


**ns**

*Depending upon concentration*

200 = 9 minutes  
400 = 16 minutes  
800 = 30 minutes  
1800 = 1 hour  
4000 = 2.2 hours  
10000 = 5.5 hours  
18000 = 10.0 hours

Each scan is 2 seconds.

Calculate total experiment time [expt



**expt**

calculates total experiment time.

**ii**

keyboard entry

May hear a **CLICK** as Probe switches in selecting nucleus slide-bar tuning positions. **Wait for II to finish.**

**zg**


keyboard entry

## Sensitivity Comparison:

DPX-300 = Av-300

DRX-400 can achieve the same Signal to Noise in half (1 / 2) the time.

**2 hours on 300 = 1 hour on 400**



# $^{13}\text{C}$ Acquisition Continued.....

OPTIONAL  
Commands  
During  $^{13}\text{C}$   
Acquisition.

**tr**  
(Optional)

The "tr" command, transfers the data from the internal console memory to the PC so that the Software can access the data while the Acquisition continues to collect additional scans.


**efp**  
(Optional)

If the Signal-to-Noise is good where you can see all your carbon peaks.....  
you may "halt" the acquisition or let it **continue** to finish collecting **ns** scans

**halt**  
(Optional)

**efp**

*Phase as needed, calibrate reference peak, etc*



# $^{13}\text{C}$ DEPT-135 Acquisition Continued.....

File....  
 ...new  
 EXPNO: 3  
 Experiment:  
 AA\_\_13C\_\_DEPT135


Prepare for a new experiment by creating a new data set and initializing its NMR parameters according to the selected experiment type.

NAME	js-072505.meod.400
EXPNO	3
PROCNO	1
DIR	D:\topspin
USER	snyder
Solvent	MeOD
Experiment	AA_13C_DEPT135
TITLE	Commercial bought sample DEPT-25-July-2005 Snyder

OK Cancel

EXPNO = 3


ns




Calculate total experiment time [expt]

Each scan is 5.6 seconds.  
 DEPT-135 setup ns such that total acquisition time is 1/2 of time needed for the carbon spectra.  
 If you ran the carbon spectra for 2 hours, then set up the DEPT-135 to run for 1 hour.

zg



efp




# $^{13}\text{C}$ & DEPT-135 run sequentially at setup time

**OPTIONAL way to Acquire the  $^{13}\text{C}$  and DEPT-135 experiments**  
**OPTIONAL way to Acquire the  $^{13}\text{C}$  and DEPT-135 experiments**  
**OPTIONAL way to Acquire the  $^{13}\text{C}$  and DEPT-135 experiments**

Create expno 2  $^{13}\text{C}$  parameters                      set ns            check acquisition time  
Create expno 3 DEPT-135 parameters                set ns            check acquisition time

Display expno 2

Click on  button

Enter in **2** for number of experiments to run.

Automation program will start on the current expno. [2]

.....When 'ns' scans are completed, or 'halt' command is given,

.....TopSpin will then read in the next incremented expno [3]

.....and start (zg) the next acquisition.

**OPTIONAL way to Acquire the  $^{13}\text{C}$  and DEPT-135 experiments**  
**OPTIONAL way to Acquire the  $^{13}\text{C}$  and DEPT-135 experiments**  
**OPTIONAL way to Acquire the  $^{13}\text{C}$  and DEPT-135 experiments**