

# <sup>1</sup>H-<sup>1</sup>H COSY Setup Bruker DPX-300 215 & Av-300 83 & DRX-400 83


on the HP Windows-XP pc computers {83 & 215} Bruker TopSpin © vers. 1.3



{This handout does not utilize Gradients, so can be run on any Bruker Avance System}


##### Handout is Still under Construction ##### Handout is Still under Construction #####

- Verify that previous acquisitions have finished [halt]
- Change samples
- Lock [lock]
- Browse to read in 1H file [browse]
- Create new datasetname [new]
- DRX-400:** Tune the <sup>1</sup>H coil of the Probe [wobb]
  
- Shim..... with sample Spinning
- ii
- rga
- zg acquire a one-scan <sup>1</sup>H



Refer to Handout “Acquire <sup>1</sup>H in preparation for 2D in TopSpin” for detailed steps for instructions below.

- Expand to desired spectral ‘Observe’ range
- Click on ‘Grid’ button for ease of locating the Spectrum’s center  
{we do not want O1 center to be on top of a peak}. With the shift-sideways button, avoid having peak at the center point }



- Click on “SW-SFO1” button to define new optimized Spectral Window for 2D.....{Set TD to 32k}



td 32k

zg

- Acquire a good S/N <sup>1</sup>H Spectra to use for Projection along 2D axis.

- Turn OFF Spinning
- Create new EXPNO 2

Spin = OFF

- Do a <sup>1</sup>H Pulse Calibration



- In EXPNO 2
- NS = 1 DS = 0 PULPROG = zg
- zg efp Touch up phase if needed.
- P1 = 20
- zg;efp Peaks should be negative phase now
- Adjust P1 until peaks are at the 360-degree null ~29<sub>usec</sub> @ PL1 = -3 dB

- Input P0 = P1 = 90 degree pulse { ¼ 360-degree null }
- Input P2 = 180 degree pulse { ½ 360-degree null }
- Input D1 = 4 sec

Spin = OFF




**Spin = OFF**




set cosy params from 1D

**PULPROG = cosyqf**  
**TD = 4k** {AQ ~ = 0.3 to 0.5 sec}  
**NS = 4, 8, 16,** multiples of 16  
**DS = 16**  
 Set **RG** to 1/2 the value from **rga** in EXPNO 1

Same as **eda** display in XwinNMR




**PULPROG = cosyqf**  
**FnMODE = QF**



**PULPROG = cosyqf**  
**FnMODE = QF**  
**TD {F2} = 4096** or 2048  
**TD {F1} = 256** {maybe 300, 360, 400, 512}  
**NS = 2, 4, 8, 16,**  
**DS = 16**  
  
**ND\_010 = 1**  
**IN\_010 = 2 x DW** usec



Scroll-down to check.... **NUC1 = 1H = 1H**




# <sup>1</sup>H-<sup>1</sup>H COSY Setup Setup Continued.....



Same as **ased** display in XwinNMR

**COSY:**  
 16 scans x 256 increments = 5 hrs 14 min.  
 8 scans x 256 increments = 2 hrs 37 min.  
 4 scans x 256 increments = 1 hrs 20 min.


Check **expt** time



**si**



**zg**



COSY FID

