

TopSpin 1.0 © Bruker BioSpin

This is the NMR Acquisition / Processing Software that will replace XwinNMR



Double click on TopSpin icon to start up software.

Your data will need to be in the directory

C:\Bruker\TOPSPIN\data\{groupname}\nmr

```
TOPSPIN 1.0
Setting environment for using Microsoft Visual C++ tools.
waiting for CPR server ...
Setting environment for using Microsoft Visual C++ tools.
CPR : Path to prog : "C:\Bruker\TOPSPIN\prog"
CPR : Path to exp : "C:\Bruker\TOPSPIN\exp"
CPR : Path to conf : "C:\Bruker\TOPSPIN\conf"
CPR : waiting for PLEXlm license
CPR : Your PLEXlm license is valid until 31-aug-2003
CPR server is up.
Starting java ...
```

