

# Create New <sup>1</sup>H filename in TopSpin



Some filename suggestions....

Do not use the same-exact filename on more than one spectrometer. As you save data, want to avoid accidentally overwriting data.

Filename examples.....Could include date info, research-log info, experiment run, spectrometer info, or reaction info etc.

js.071505  
js.071505-1  
tpp.fr3.cdcl3-  
log214.expt7.400  
log214.expt7.bc  
log214.expt7.cono  
log214.expt7.inv  
log214.expt7.600  
phen\_br\_150C\_36hrs

Do not use "/" or "\" in filename!

File....  
...new

or...




**Cntrl + N**

## 1. Change Filename [NAME]

## 2. Change Title

**Tip:** May want to include....  
what sample is  
solvent info  
spectrometer  
date

temperature  
pH  
Nucleus  
etc...



## OPTIONS:

As needed.....

Can update Solvent

Can read in Parameter file  
[Experiment]

